
# TIGR® Matrix Surgical Mesh

Nätimplantat / Chirurgisches Netz / Malha Cirúrgica / Malla Quirúrgica / Chirurgická Síťka  
Maglia Chirurgica / Chirurgisch Gaas / Kirurgisk Net / Maille Chirurgicale  
Siátka Chirurgiczna / Kirurgisk Nett/ Cerrahi Mesh / Kirurginen verkko / Χειρουργικό πλέγμα  
Kirurška mrežica / Hirurška mrežica / شبكة جراحية / Хирургическая мрежа

## Instructions for use

Bruksanvisning / Gebrauchsanweisung / Instruções de Uso / Instrucciones de Uso / Návod k Použití  
Istruzioni per l'uso / Instructies voor het Gebruik / Brugsanvisning / Mode d'emploi  
Instrukcja Użytkowania / Instruksjoner for Bruk / Kullanım Talimatları / Käyttöohjeet / Οδηγίες Χρήσης  
Upute za uporabu / Uputstvo za upotrebu / تعليمات للاستخدام / Инструкции за употреба

- EN** PAGE 3 FOR ENGLISH  
SEE INSIDE OF LEAFLET FOR SYMBOL EXPLANATION IN ENGLISH
- SV** SIDA 4 FÖR SVENSKA  
SE INSIDAN AV BOKEN FÖR SYMBOLFÖRKLARING PÅ SVENSKA
- DE** SEITE 5 FÜR DEUTSCH  
SIEHE PACKUNGSBEILAGE FÜR EINE ERKLÄRUNG DER SYMBOLE AUF DEUTSCH
- PT** PÁGINA 6 PARA PORTUGUÊS  
VEJA PARA DENTRO PARA A EXPLICAÇÃO DOS SÍMBOLOS EM PORTUGUÊS
- ES** PÁGINA 7 PARA ESPAÑOL  
VER EL INTERIOR DEL FOLLETO PARA LA EXPLICACIÓN DE LOS SÍMBOLOS EN ESPAÑOL
- CS** STRANA 8 V ČESKÉ  
VČESKÉ VYSVĚTLIVKY K SYMBOŮM JSOU UVEDENY V PŘÍBALOVÉ INFORMACI
- IT** PAGINA 9 PER L'ITALIANO  
PER LA SPIEGAZIONE DEI SIMBOLI IN ITALIANO, SI VEDA IL FOGLIETTO ILLUSTRATIVO
- NL** PAGINA 10 VOOR NEDERLANDS  
ZIE DE BINNENKANT VAN DE BIJSLUITER VOOR NEDERLANDSE TOELICHTING VAN DE SYMBOLEN
- DA** SIDE 11 I DANSK  
SE INDERSIDEN AF INDLÆGSSEDLEN FOR SYMBOLFORKLARING I DANSK
- FR** PAGE 12 POUR LE FRANÇAIS  
VOIR L'INTÉRIEUR DU FEUILLET POUR CONNAÎTRE L'EXPLICATION DES SYMBOLES EN FRANÇAIS
- PL** STRONA 13 DLA POLSKI  
ZOBACZ WNĘTRZE ULOTKI DLA SYMBOLI OBJAŚNIENYCH W JĘZYKU POLSKIM
- NO** SIDE 14 FOR NORSK  
SE INNSIDEN AV PAKKEN FOR FORKLARING AV SYMBOLER PÅ NORSK
- TR** SAYFA 15 TÜRKÇE  
SEMBOŖLÜN TÜRK AÇIKLAMASI İÇİN PROSPEKTÜSE BAKIN
- FI** SIVU 16 SUOMEKSI  
SUOMEKSI KATSO PAKKAUSSELOSTEEN SISÄSIVUILLA SYMBOLIEN SELITYKSET SUOMEKSI
- EL** ΣΕΛΙΔΑ 17 ΓΙΑ ΕΛΛΗΝΙΚΑ  
ΔΕΙΤΕ ΤΟ ΕΣΩΤΕΡΙΚΟ ΤΟΥ ΦΥΛΛΑΔΙΟΥ ΓΙΑ ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ ΣΤΑ ΕΛΛΗΝΙΚΑ
- SR** STRANA 18 ZA SRPSKI JEZIK  
OBJAŠNJENJE SIMBOLA NA SRPSKOM JEZIKU POTRAŽITE U UPUTSTVU ZA UPOTREBU
- HR** STRANA 19 ZA HRVATSKI  
ZA OBJAŠNJENJE SIMBOLA POGLEDATI UPUTSTVA ZA UPORABU
- AR** الصفحة 20 للغة العربية  
اطلع على النشرة لتوضيح الرموز باللغة العربية
- BG** СТРАНИЦА 21 ЗА БЪЛГАРСКИ  
ВИЖТЕ ВЪТРЕШНАТА СТРАНА НА ЛИСТОВКАТА ЗА ОБЯСНЕНИЕ НА СИМВОЛИТЕ НА АНГЛИЙСКИ ЕЗИК

**DEVICE DESCRIPTION**

TIGR® Matrix Surgical Mesh is knitted from two different synthetic resorbable fibers, possessing different degradation characteristics. The fast-resorbing fiber, making up approximately 40% of the matrix by weight, is a copolymer of glycolide, lactide, and trimethylene carbonate. The slow-resorbing fiber, making up approximately 60% of the matrix by weight, is a copolymer of lactide, and trimethylene carbonate. Both fibers degrade by bulk hydrolysis once implanted, resulting in a decreasing strength retention followed by mass loss of the fibers. In vitro testing showed that the fast-resorbing fiber (glycolide, lactide, and trimethylene carbonate) loses its mechanical strength after 2 weeks and in vivo studies in the abdominal wall of sheep showed that the fast-resorbing fiber is fully absorbed after 4 months. The same in vitro testing showed that the slow-resorbing fiber (lactide, and trimethylene carbonate) maintains its mechanical strength for 6 months and in vivo studies in the abdominal wall of sheep indicated that the slow-resorbing fiber is absorbed after approximately 36 months.

**INDICATIONS FOR USE**

TIGR® Matrix Surgical Mesh is intended for use in reinforcement of soft tissue where weakness exists in procedures involving repair of indirect inguinal hernias and ventral hernias, in prophylactic use to reinforce the midline sutures, and in reconstructive breast surgery for both prepectoral and submuscular surgical procedures.

**CONTRAINDICATIONS**

Not suitable for reconstruction of cardiovascular defects. TIGR® Matrix Surgical Mesh must always be separated from the abdominal cavity by peritoneum. Not suitable for the repair of direct inguinal hernias.

**WARNINGS**

1. Do not use if the outer or inner package has been damaged and/or any sterile or moisture barrier is not intact.
2. Do not use after the expiration date – the biodegradable components may not perform adequately.
3. The use of any synthetic mesh or patch in a contaminated or infected wound could lead to fistula formation and/or extrusion of the mesh and it is not recommended.
4. For single use only. Do not use if outer or inner package has been opened prior to initial intended use. Discard all unused portions of the device. Do not sterilize. The device might not perform adequately due to degradation.
5. The safety and effectiveness of TIGR® Matrix Surgical Mesh have not been established for use with resorbable fixation devices.
6. The safety and effectiveness of TIGR® Matrix Surgical Mesh have not been established for uroynecological use. Refer to safety communications from the FDA and from UK's National Institute for Health and Clinical Excellence (NICE) for guidance.
7. The safety and effectiveness of TIGR® Matrix Surgical Mesh have not been established for use in tendon repair.
8. Because TIGR® Matrix Surgical Mesh is fully resorbable, it should not be used in repairs where permanent support from the mesh is required.

**PRECAUTIONS**

1. This device is restricted for use by or on the order of a physician.
2. Carefully check that the packaging is undamaged and unopened and that the sterile barrier is intact before use.
3. The mesh should be large enough to extend beyond the margin of the defect.
4. Infections should be treated according to acceptable surgical practice to minimize the need for removal of the mesh.

**ADVERSE REACTIONS**

Possible adverse reactions with the mesh are those typically associated with any implantable prosthesis, including, but not limited to, infection, inflammation, extrusion, erosion, adhesion, fistula formation, seroma formation, hematoma, and recurrence of the hernia or tissue defect. Known adverse reactions following breast reconstructive surgery with or without mesh includes, but not limited to, implant displacement, infection, inflammation, hematoma, seroma formation, necrosis, capsular contraction and implant rupture and in rare cases breast Implant associated anaplastic large cell lymphoma (BIA-ALCL).

**CLINICAL BENEFITS**

Depending on the indication for use, the following clinical benefits can be expected:

- Freedom from hernia-related symptoms.
- Prevention of surgery-related incisional hernia.
- Short- and long-term successful subpectoral and prepectoral breast reconstruction.

The magnitude of the clinical benefits related to TIGR® Matrix Surgical Mesh is comparable with that of other commercially available surgical meshes.

**PREPARATION FOR USE**

1. Open the outer aluminum foil and remove the inner pouch containing the product. Be aware that the inner pouch is non-sterile on the outside.
2. Carefully open the inner pouch containing the mesh.
3. Aseptically remove the mesh from the inner pouch using sterile, gloved hands or sterile forceps, and place the mesh in the sterile field.

**DIRECTIONS FOR USE**

1. Prepare the implantation site using standard surgical techniques.
2. Trim TIGR® Matrix Surgical Mesh to allow an adequate overlap of the defect area or to fit the size of the implant when used for breast reconstruction.
3. Implant TIGR® Matrix Surgical Mesh according to currently accepted surgical mesh procedures. Device may be used in a dry state.
4. Fixate TIGR® Matrix Surgical Mesh with sutures according to currently accepted surgical practices.
5. Affix the traceability label in the patient's medical record and provide the implant card to the patient.
6. If a serious incident should occur in relation to the device it shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

**STORAGE, PACKAGING AND DISPOSAL**

1. Store in a cold dry place away from moisture and direct heat inside the product box with unopened and undamaged aluminum pouch, instruction for use and implant card.
2. Sterile in unopened and undamaged package with sterile barrier intact.
3. A traceability label which identifies the lot number of the product is enclosed in every package for placement in the patient's medical record.
4. An implant card which identifies the product is enclosed in every package.
5. Dispose of contaminated units, components, and packaging materials in accordance with standard hospital procedures, universal precautions for biohazardous waste, and applicable local, state, and federal laws.

**DISCLAIMER OF WARRANTY**

Although TIGR® Matrix Surgical Mesh (hereinafter referred to as "product") has been manufactured under carefully controlled conditions; Novus Scientific AB (hereinafter called Novus) has no control over the conditions under which the product is used. Novus, therefore, disclaims all warranties, both expressed and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Novus shall not be liable to any person or entity for any medical expenses or any direct, indirect, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise.

No person has any authority to bind Novus to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law, including the Federal Drug, and Cosmetic Act. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

**EXPLANATION OF SYMBOLS**

1. Caution:  
This device is restricted for use by or on the order of a physician.
2. Consult Instructions for Use.
3. Use before date.
4. Sterilized using ethylene oxide.
5. Single sterile barrier system with protective packaging outside. (Used on the aluminum pouch.)
6. For single use only.
7. Catalogue number.
8. Lot number.
9. Date and country of manufacture.
10. Quantity.
11. Storage conditions.
12. Manufacturer.
13. Complies with the Medical Device Regulation 2017/745.
14. Do not sterilize.
15. Do not use if package is damaged.
16. Medical device.
17. Unique device identifier.
18. Patient information website.
19. Person identification.
20. Date.
21. Health care center or doctor.
22. Single sterile barrier. (Used on the Tyvek pouch.)

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TIGR® is a registered trademark for Novus Scientific. The product is protected by EP1674048; EP1870056; EP2002800; US 8,016,841; US 8,083,755; US 8,313,499; US 8,906,047; US 9,566,370; US 9,668,847; US 9,717,825; US 9,750,854; US 10,342,653; AU2008202439; CN 10140748 & JP5203808.