
TIGR[®] Matrix Surgical Mesh

Nätimplantat / Chirurgisches Netz / Malha Cirúrgica / Malla Quirúrgica / Chirurgická Síťka
Maglia Chirurgica / Chirurgisch Gaas / Kirurgisk Net / Maille Chirurgicale
Siatka Chirurgiczna / Kirurgisk Nett/ Cerahi Mesh / Kirurginen verkko / Χειρουργικό πλέγμα
Kirurška mrežica

Instructions for use

Bruksanvisning / Gebrauchsanweisung / Instruções de Uso / Instrucciones de Uso / Návod k Použití
Istruzioni per l'uso / Instructies voor het Gebruik / Brugsanvisning / Mode d'emploi
Instrukcja Użytkowania / Instruksjoner for Bruk / Kullanım Talimatları / Käyttöohjeet / Οδηγίες Χρήσης
Upute za uporabu

- (EN) PAGE 2 FOR ENGLISH
SEE INSIDE OF LEAFLET FOR SYMBOL EXPLANATION IN ENGLISH
- (SV) SIDA 3 FÖR SVENSKA
SE INSIDAN AV BOKEN FÖR SYMBOLFÖRKLARING PÅ SVENSKA
- (DE) SEITE 4 FÜR DEUTSCH
SIEHE PACKUNGSBEILAGE FÜR EINE ERKLÄRUNG DER SYMBOLE AUF DEUTSCH
- (PT) PÁGINA 5 PARA PORTUGUÊS
VEJA PARA DENTRO PARA A EXPLICAÇÃO DOS SÍMBOLOS EM PORTUGUÊS
- (ES) PÁGINA 6 PARA ESPAÑOL
VER EL INTERIOR DEL FOLLETO PARA LA EXPLICACIÓN DE LOS SIMBOLOS EN ESPAÑOL
- (CS) STRANA 7 V ČESKÉ
ČESKÉ VYSVĚTLIVKY K SYMBOLŮM JSOU UVEDENY V PŘÍBALOVÉ INFORMACI
- (IT) PAGINA 8 PER L'ITALIANO
PER LA SPIEGAZIONE DEI SIMBOLI IN ITALIANO, SI VEDA IL FOGLIETTO ILLUSTRATIVO
- (NL) PAGINA 9 VOOR NEDERLANDS
ZIE DE BINNENKANT VAN DE BIJSLUITER VOOR NEDERLANDSE TOELICHTING VAN DE SYMBOLEN
- (DA) SIDE 10 I DANSK
SE INDERSIDEN AF INDLÆGSSEDLEREN FOR SYMBOLFORKLARING I DANSK
- (FR) PAGE 11 POUR LE FRANÇAIS
VOIR L'INTÉRIEUR DU FEUILLET POUR CONNAÎTRE L'EXPLICATION DES SYMBOLES EN FRANÇAIS
- (PL) STRONA 12 DLA POLSKICH
ZOBACZ WNĘTRZE ULOTKI DLA SYMBOLI OBJAŚNIONYCH W JĘZYKU POLSKIM
- (NO) SIDE13 FOR NORSK
SE INNSIDEN AV PAKKEN FOR FORKLARING AV SYMBOLER PÅ NORSK
- (TR) SAYFA14 TÜRKÇE
SEMBOLOUN TÜRK AÇIKLAMASI İÇİN PROSEKTÜSE BAKIN
- (FI) SIVU 15 SUOMEKSI
SUOMEKSI KATSO PAKKAUSSELOSTEEN SISÄSIVUILTA SYMBOLIEN SELITYKSET SUOMEKSI
- (EL) ΣΕΛΙΔΑ16 ΓΙΑ ΕΛΛΗΝΙΚΑ
ΔΕΙΤΕ ΤΟ ΕΣΩΤΕΡΙΚΟ ΤΟΥ ΦΥΛΛΑΔΙΟΥ ΓΙΑ ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ ΣΤΑ ΕΛΛΗΝΙΚΑ
- (HR) STRANA 18 ZA HRVATSKI
ZA OBJAŠNJENJE SIMBOLA POGLEDATI UPUTSTVA ZA UPORABU

Instructions for Use – ENGLISH

EN

DEVICE DESCRIPTION

TIGR® Matrix Surgical Mesh is knitted from two different synthetic resorbable fibers, possessing different degradation characteristics. The fast-resorbing fiber, making up approximately 40% of the matrix by weight, is a copolymer of glycolide, lactide, and trimethylene carbonate. The slow-resorbing fiber, making up approximately 60% of the matrix by weight, is a copolymer of lactide, and trimethylene carbonate. Both fibers degrade by bulk hydrolysis once implanted, resulting in a decreasing strength retention followed by mass loss of the fibers. In vitro testing showed that the fast-resorbing fiber (glycolide, lactide, and trimethylene carbonate) loses its mechanical strength after 2 weeks and in vivo studies in the abdominal wall of sheep showed that the fast-resorbing fiber is fully absorbed after 4 months. The same in vitro testing showed that the slow-resorbing fiber (lactide, and trimethylene carbonate) maintains its mechanical strength for 6 months and in vivo studies in the abdominal wall of sheep indicated that the slow-resorbing fiber is absorbed after approximately 36 months.

INDICATIONS FOR USE

TIGR® Matrix Surgical Mesh is intended for use in reinforcement of soft tissue where weakness exists, in procedures involving the repair of hernias and abdominal wall defects, abdominal wall reinforcement and muscle flap reinforcement.

CONTRAINDICATIONS

Not suitable for reconstruction of cardiovascular defects. TIGR® Matrix Surgical Mesh must always be separated from the abdominal cavity by peritoneum. Not for use following planned intra-operative or accidental opening of the gastrointestinal tract. Use in these cases may result in contamination of the mesh, which may lead to infection. Not suitable for the repair of direct inguinal hernias.

WARNINGS

1. Do not use if the outer or inner package has been damaged and/or any sterile or moisture barrier is not intact.
2. Do not use after the expiration date – the biodegradable components may not perform adequately.
3. Do not use on contaminated and/or infected wounds.
4. For single use only. Do not use if outer or inner package has been opened prior to initial intended use. Discard all unused portions of the device. Do not resterilize. The device might not perform adequately due to degradation.
5. The safety and effectiveness of TIGR® Matrix Surgical Mesh have not been established for use with resorbable fixation devices.
6. The safety and effectiveness of TIGR® Matrix Surgical Mesh have not been established for urogynecological use. Refer to safety communications from the FDA and from UK's National Institute for Health and Clinical Excellence (NICE) for guidance.
7. The safety and effectiveness of TIGR® Matrix Surgical Mesh have not been established for use in tendon repair.
8. Because TIGR® Matrix Surgical Mesh is fully resorbable, it should not be used in repairs where permanent support from the mesh is required.

PRECAUTIONS

1. This device is restricted for use by or on the order of a physician.
2. Carefully check that the packaging is undamaged and unopened and that the sterile barrier is intact before use.
3. The mesh should be large enough to extend beyond the margin of the defect.
4. Infections should be treated according to acceptable surgical practice to minimize the need for removal of the mesh.

ADVERSE REACTIONS

Possible adverse reactions with the mesh are those typically associated with any implantable prosthesis, including, but not limited to, infection, inflammation, extrusion, erosion, adhesion, fistula formation, seroma formation, hematoma, and recurrence of the hernia or tissue defect.

PREPARATION FOR USE

1. Open the outer aluminum foil and remove the inner pouch containing the product. Be aware that the inner pouch is non-sterile on the outside.
2. Carefully open the inner pouch containing the mesh.
3. Aseptically remove the mesh from the inner pouch using sterile, gloved hands or sterile forceps, and place the mesh in the sterile field.

DIRECTIONS FOR USE

1. Prepare the implantation site using standard surgical techniques.
2. Trim TIGR® Matrix Surgical Mesh so as to allow an adequate overlap of the defect area.
3. Implant TIGR® Matrix Surgical Mesh according to currently accepted surgical mesh procedures either open (e.g. acc to Lichtenstein, TIPP) or laparoscopic (e.g. TAPP, TEP).
4. Fixate TIGR® Matrix Surgical Mesh with sutures or staples according to currently accepted surgical practices.
5. Affix the traceability label in the patient's medical record.
6. Device may be used in dry state.

STORAGE, PACKAGING AND DISPOSAL

1. Store in a cool dry place away from moisture and direct heat.
2. Sterile in unopened and undamaged package with sterile barrier intact.
3. A traceability label which identifies the lot number of the prosthesis is enclosed in every package for placement in the patient's medical record.
4. Dispose of contaminated units, components, and packaging materials in accordance with standard hospital procedures, universal precautions for biohazardous waste, and applicable local, state, and federal laws.

DISCLAIMER OF WARRANTY

Although TIGR® Matrix Surgical Mesh (hereinafter referred to as "product") has been manufactured under carefully controlled conditions; Novus Scientific Pte. Ltd (hereinafter called Novus) has no control over the conditions under which the product is used. Novus, therefore, disclaims all warranties, both expressed and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Novus shall not be liable to any person or entity for any medical expenses or any direct, indirect, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Novus to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law, including the Federal Drug, and Cosmetic Act. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

EXPLANATION OF SYMBOLS

1. Caution: This device is restricted for use by or on the order of a physician.
2. Consult Instructions for Use.
3. Use before date.
4. Sterilized using ethylene oxide.
5. Tyvek pouch non-sterile on outside!
6. For single use only.
7. Catalogue number.
8. Lot number.
9. Date of manufacture.
10. Quantity.
11. Storage conditions.
12. Manufacturer.
13. Complies with the Medical Device Directive 93/42/EEC amended through Directive 2007/47/EC.
14. Do not resterilize.
15. Do not use if package is damaged.

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TIGR[®] is a registered trademark of Novus Scientific. The product is protected by:
Novus Scientific AB